




<b>Product</b> <b>Item number</b>	CareSystem backrest, height and sideways adjustable R6533000
<b>Colour</b>  <b>Product description</b>	<div>  </div> <p>White</p> <p><b>Discontinued 31 Dec 2015.</b> Back rest for shower chair for toilet without cistern. 150 mm height adjustable, for mounting on horizontal track.  (recommended for toilet pans with a total distance of 700 mm from wall to front)</p> <p>Back rest for shower chair for toilet without cistern. 150 mm height adjustable, for mounting on horizontal track.  (recommended for toilet pans with a total distance of 700 mm from wall to front)</p>
<b>Adjustability</b> <b>Product capacity</b>	Height and sideways adjustable Max. vertical load: 100 kg.
<b>Materials</b>  <b>Ergonomy</b> <b>Features</b>	<p>Basic material: 20 micrometre anodized aluminium</p> <p>Backrest: Polyurethane foam with marine ply core</p> <p>Backrest surface: Gum lac</p> <p>Lever: Stainless steel w/nylon knob</p> <p>Front plate/end covers, track: Polystyrene plastic (PS)</p> <p>Soft and slip-resistant material with a warm feel.</p> <p>For toilet without cistern.</p> <p>Recommended when the toilet's projection = 700 mm.</p> <p>Can be pushed along the track.</p> <p>Height adjustable 150 mm.</p>



# PRESSALIT®

<b>Product</b>	CareSystem backrest, height and sideways adjustable
<b>Item number</b>	R6533000
<b>Cleaning</b>	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

## PRESSALIT®