



# PRESSALIT®

|                                      |   |
|--------------------------------------|---|
| <b>Product</b><br><b>Item number</b> | Changing table, 800 x 1800 mm, electrically height adjustable<br>R8743000   |
| <b>Colour</b>                        | White   |
| <b>Product description</b>           | <p><b>Discontinued per 31 Dec 2015 and replaced by R8643</b></p> <p>Freestanding changing table, 800 x 1800 mm, with electric motor, with remote control, incl. transformer. Height adjustable from 650 mm - 1000 mm.<br/>Excl. mounting kit.</p> <p><b>Discontinued per 31 Dec 2015 and replaced by R8643</b></p> <p>Freestanding changing table, 800 x 1800 mm, with electric motor, with remote control, incl. transformer. Height adjustable from 650 mm - 1000 mm.<br/>Excl. mounting kit.</p> |
| <b>Adjustability</b>                 | Height adjustable   |
| <b>Product capacity</b>              | Max load: 150 kg.   |
| <b>Speed upwards</b>                 | Speed up mm: 14 mm/s<br>Speed down mm: 14 mm/s  |
| <b>Noise level</b>                   | 48 dB   |
| <b>Materials</b>                     | Lifting unit and baskets: Powder-laquered steel.<br>Worktop: Melamine   |
| <b>Ergonomy</b>                      | The smooth height adjustment ensures that it is always possible to set the table at an ergonomically correct height.  |
| <b>Features</b>                      | Easy-to-clean surfaces ensure optimum hygiene.  |
| <b>Mains connection (V)</b>          | 230   |
| <b>Actuator (motor)</b>              | 2 x LA31 Deskline (Linak), 1500N, IPX1  |
| <b>Control box</b>                   | CBD6 (Linak), IP22  |



# PRESSALIT®

|   |  |
|---|--|
| <b>Product</b>                          | Changing table, 800 x 1800 mm, electrically height adjustable  |
| <b>Item number</b>                      | R8743000   |
| <b>Wired remote control</b>             | HB3 (Linak), IP66  |
| <b>Power consumption (maximum load)</b> | Power consumption by maximum load, watt: 96 W<br>Power consumption by standby: 0,1 W   |
| <b>Cleaning</b>                         | The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®